ROYAL®
COLLECTION

čeština	1
slovenčina	11
english	21
deutsch	31
русский	41
magyar	51

čeština

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Lady Dreams®

English

Congratulations! You have made a great choice by selecting and purchasing this sex toy.

Microscopium Ben Wa balls feature an integrated innovative “trampoline effect”.

In this way subtle, “kinetic vibrations” are triggered which stimulate your pelvic floor muscles to contract and relax while also providing you with stimulation. The added value lies in that this allows you to strengthen the pelvic floor with minimum effort.

Microscopium Ben Wa balls feature a special pull-out loop, which is worn discreetly - hidden inside the body. This makes them some of the most discreet beads available on the market. So you can wear them everywhere, even when swimming, in the sauna, etc. Very easy to implement and also easily removable, thanks to the innovative loop.

Advice for beginners – in the beginning it’s quite possible that the balls may slip out; in this case treat the balls hygienically and reinsert.

A perfect tool for training the pelvic floor and improving elasticity of the vaginal muscles. Unbeatable tool used to prevent (as well as a form of treatment for) incontinence - accidental leakage of urine.

Instructions for maintaining Microscopium Ben Wa balls

1. Maintenance of the Microscopium Ben Wa balls is simple. Before and after each use, rinse the Microscopium

balls and disinfect them with a disinfectant designed specifically for treating sex toys. This ensures you protect your health while extending the lifespan of your sex toy. Never use petrol or alcohol for cleaning as they may damage both the surface and the product.

2. In order to achieve the most pleasant feelings and the maximum satisfaction, apply the sufficient amount of a lubricant to the surface of the Microscopium Ben Wa balls. Only use oil-free lubricants, which do not affect the material.

Instructions for using Microscopium Ben Wa balls

1. Choose the most appropriate position for inserting the Microscopium balls - for beginners the most popular position for insertion is when lying down or with one leg raised. Lubricate the balls before introduction – then the actual introduction will be easy and pleasant.

2. As with tampons, you should soon get used to this.

3. Wear the Microscopium balls for at least 15 minutes; for improved exercise, wearing for several hours is ideal; however do not exceed a maximum of 8 hours.

4. You will also feel some stimulation with the gradual

„up and down“ application, i.e. by pulling and repeated pushing in the Microscopium balls in regular cycles.

5. **IMPORTANT:** The Microscopium balls can never disappear in the body – simply put, there is no place for them to disappear to.

6. If you feel pain even when introducing very gently - consult your doctor.

7. Keep your Microscopium balls in the clean, cool and dry condition. Protect the device from direct sunlight!

8. Store your Microscopium balls separately rather than together with other products of dissimilar materials in order to avoid chemical reactions.

Exercising tips

You can easily train, for example, while performing physical activities - such as walking, aerobics, rollerblading, horseback riding and dancing etc. For comfortable training the activity can be replaced by rhythmic movements such as jumping on a ball or rocking chair.

Maximum contractions

After introducing the Microscopium balls just lie on your back and bend your knees. Contract your pelvic floor

muscles, as if you wanted to hold in urine and intestinal gases, pull in your belly and slowly exhale through your mouth. Hold the contraction for 5 seconds and then add three quick contractions in turn. Relax the pelvic floor and abdomen: Breathing will occur automatically. Do three series of ten very strong contractions of the pelvic floor muscles. Have ten-second rests between individual contractions and one-minute rests between the individual series. Relax all muscles which could replace pelvic floor muscle activity (gluteal muscles and inner thigh muscles) as much as you can. Practice once or twice a day for 5 days a week. To make sure that you are contracting the proper muscles, always introduce the Microscopium balls before exercising. If you feel that the Microscopium balls are being pulled inside the vagina, this means that the pelvic floor muscle contractions are being carried out properly. If, however, the Microscopium balls are being ejected from the vagina on contraction, this means that instead of muscle contractions, you are pressing (as when defecating). Attention! It is very important to properly understand the movement before repeating. Pushing out (as when defecating) might weaken the pelvic floor even further. If you fail to do the pelvic floor muscles, try the following exercise.

Wave

After introducing the Microscopium balls just lie on your

back and bend your knees. Pinch the rectum, as if you wanted to stop intestinal gases, while exhaling through your mouth and retracting your abdomen. Try contracting your pelvic floor muscles from the rectum toward the vagina. When the rectum and vagina are contracted, maintain the contraction for 5 seconds and then relax the pelvic floor and the abdomen. Repeat several times and check the direction, in which the Microscopium balls draw. Once you have mastered this exercise, try maximum contraction of the pelvic floor muscles once again.

Note: If you really want to strengthen your pelvic floor muscles, the contraction must be as strong as possible, which means that you have to do it in a concentrated and intense manner.

Locking the dam

The exercise consists of contracting your pelvic floor prior to any physical activity which increases pressure on the dam, such as coughing, sneezing, or lifting a child, and also during this activity. Get used to performing dam locking during any physical exertion.

Geisha

During sexual intercourse when being penetrated, contract the pelvic floor muscles and squeeze the vagina. Hold the contraction for a few seconds and then release

the grip. Repeat twice to three times in a row. Your partner will let you know whether your pelvic floor muscles are stronger.

Targeted contraction exercises of the pelvic floor
- active strengthening

Introduce the Microscopium Ben Wa balls according to the instructions above. Get into a slight crouch and pull the pelvic floor muscles. Then hold the loop and try to pull the ball out, while avoiding pulling it through active muscle contraction. Start training for a shorter period and then steadily built it up. Practice in this way on a daily basis at most for 15 minutes.

Exercise warning

Even vaginal muscles exercises may cause pain and muscle fatigue. Therefore, make sure that you do not exercise too long in the beginning. It is better to gradually build up your exercise periods.

Warnings and general instructions

1. Keep out of reach of children. Store in a clean, cool and dry environment. Protect from direct sunlight and do not expose to extreme temperatures.
2. Do not use your Microscopium Ben Wa balls if you have diseased and irritated skin and do not use it on

infected parts of your skin. The Microscopium Ben Wa balls are not designed and intended as a medical aid.

3. Read the instructions on the packaging regarding the material to avoid allergic reactions (for example latex).

4. If you feel any awkward tension or irritation of the skin while using the Microscopium balls, discontinue use immediately.

5. Do not bend or stretch the Microscopium Ben Wa balls too strongly. Do not apply force.

6. If the Microscopium Ben Wa balls are damaged do not use them any more.

7. The supplier does not assume any responsibility for the product if it is not used properly.

Environmental protection warning

Do not throw packaging material and used devices away; submit them for recycling instead.

Warranty

We offer a two-year warranty for Microscopium Ben Wa balls starting from the date of their original purchase. The warranty does not apply if destroyed, misused

or due to any other external cause. Any attempt to repair the product yourself will void the warranty. The customer is obliged to present the product purchase receipt when claiming a refund.

Thank you for your purchase. For other interesting information, discounts, competitions, news, etc. see www.eroticcity.cz.

Producer

City Realex, a.s., Tenisová 16, Praha 10,
Czech Republic

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